



Be Right™



Solución estándar de calcio, 100 mg/L como Ca (NIST), 500 mL

de producto:

4457649

USD Precio:

Contacto Hach

Disponible

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

Para calibrar los medidores y comprobar la exactitud de las determinaciones de calcio. Botella de 500 mL.

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Especificaciones

Cantidad:	500 mL
Concentración:	100 mg/L Ca
Condiciones de almacenamiento:	10 - 25 °C (protect from light)
Notas especiales:	Observed shelf life after opening, following recommended storage and use: 3 months
Parámetro:	Calcium
Plataforma:	Solución
Tipo de envase:	Bottle